

RMassBank: Run-through the Principles and Workflow in R

Emma Schymanski

Michael Stravs, Heinz Singer & Juliane Hollender

Eawag, Dübendorf, Switzerland



Steffen Neumann, Erik Müller: *IPB Halle, Germany*

Tobias Schulze: *UFZ Leipzig, Germany*

RMassBank Questions: massbank@eawag.ch

Demonstration Overview

- Install (and update) RMassBank and associated programs
- Run RMassBank for trial data
- Get a basic understanding for the workflow
- Understand the manual checking required
 - “Fail peaks” – checking these in raw data
 - Automatic annotation with CTS – editing this data
- Generate records with RMassBank
- Upload these records to MassBank



MassBank Record: EA067410

[Home](#) | [Spectrum](#) | [Query](#) | [Data](#) | [Substances](#) | [Browse](#) | [Database](#) | [Index](#) | [MassBank-ID](#)

Verapamil, LC-ESI-ITFT, MS2, 45%; R=16000; [M+H]⁺

Mass Spectrometry:

Chemical Structure:

ACQUISITION: EA067410
 RECORD_TITLE: Verapamil; LC-ESI-ITFT; MS2; 45%; R=16000; [M+H]⁺
 DATE: 2011.09.04
 AUTHOR: Steffy M, Schymanski E, Singer M, Department of Environmental Chemistry, Eawag
 LICENSE: <https://massbank.eur.de/massbank/files/license.html>
 COPYRIGHT: Copyright (C) 2011 Eawag, Dübendorf, Switzerland
 COMMENT: CONFERENCE: wuolabs congress
 COMMENT: EAWAG_UPLOAD_ID: 474

Starting RMassBank

Open up RStudio (or R console)

- o And open up “RMassBank_Demo.R”

The screenshot shows the RStudio interface with the RMassBank_Demo.R script open in the editor and the R console displaying the output of the script. The script contains the following R code:

```

1 # Interactive RMassBank Demonstration
2 # E. Schymanski, Eawag, 27/11/2012
3
4 # First, check that you have the latest stable RMassBank installed
5 source("http://bioconductor.org/biocLite.R")
6 # Identify packages requiring update
7 old.packages(repos=biocinstallRepos())
8 # update all packages that are out of date
9 biocLite(character(), ask=FALSE)
10 # or just install RMassBank (does the rest for you)
11 source("http://bioconductor.org/biocLite.R")
12 biocLite("RMassBank")
13 biocLite("RMassBankData")
14 # Careful, if you haven't updated for a while this can
15 # require update of several related packages (several minutes)
16

```

The R console output shows the following messages:

```

Natural language support but running in an English locale

R is a collaborative project with many contributors.
Type 'contributors()' for more information and
'citation()' on how to cite R or R packages in publications.

Type 'demo()' for some demos, 'help()' for on-line help, or
'help.start()' for an HTML browser interface to help.
Type 'q()' to quit R.

Loading required package: rcom
Loading required package: rscproxy
[Workspace loaded from ~/.RData]

Warning message:
Class "C++Object" is defined (with package slot 'Rcpp') but no metadata object found
to revise subclass information---not exported? Making a copy in package '.GlobalEnv'
>

```

The Packages pane on the right shows a list of installed and available packages, including abind, acepack, actuar, ada, adaptMCMC, AER, akima, AlgDesign, ALL, ALL, alr2, amap, anchors, animation, and ape.

Installing / Updating RMassBank

To check if you need to update RMassBank:

```

1 # Interactive RMassBank Demonstration
2 # E. Schymanski, Eawag, 27/11/2012
3
4 # First, check that you have the latest stable RMassBank installed
5 source("http://bioconductor.org/biocLite.R")
6 # Identify packages requiring update
7 old.packages(repos=biocinstallRepos())
8 # update all packages that are out of date
9 biocLite(character(), ask=FALSE)
10 # or just install RMassBank (does the rest for you)
11 source("http://bioconductor.org/biocLite.R")
12 biocLite("RMassBank")
13 biocLite("RMassBank")
14 # Careful, if you get an error message, it might be a problem with the
15 # require up
16 # if you get an error message, it might be a problem with the

```

Console ~ / ↻

```


> source("http://bioconductor.org/biocLite.R")
BiocInstaller version 1.2.1, ?biocLite for help
> old.packages(repos=biocinstallRepos())

```

Package	LibPath
"bitops"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"Cairo"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"coda"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"fingerprint"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"RANN"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"rgl"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"rscproxy"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"XML"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"zoo"	"C:/Users/schymaem/Documents/R/win-library/2.14"
"acepack"	"C:/Program Files/R/R-2.14.1/library"
"ape"	"C:/Program Files/R/R-2.14.1/library"
"bitops"	"C:/Program Files/R/R-2.14.1/library"

Installing / Updating RMassBank

To install RMassBank:



```
1 # Interactive RMassBank Demonstration
2 # E. Schymanski, Eawag, 27/11/2012
3
4 # First, check that you have the latest stable RMassBank installed
5 source("http://bioconductor.org/biocLite.R")
6 # Identify packages requiring update
7 old.packages > biocLite("RMassBank")
8 # update
9 biocLite("RMassBank")
10 # or just
11 source("http://bioconductor.org/biocLite.R")
12 biocLite("RMassBank")
13 biocLite("RMassBank")
14 # Careful
15 # require
16 # if you
```

BioC_mirror: 'http://www.bioconductor.org'
Using R version 2.14, BiocInstaller version 1.2.1.
Installing package(s) 'RMassBank'
Installing package(s) into 'C:/Users/schymaem/Documents/R/win-library/2.14'
(as 'lib' is unspecified)
Old packages: 'bitops', 'Cairo', 'coda', 'fingerprint', 'RANN', 'rgl', 'rscproxy',
'XML', 'zoo', 'acepack', 'ape', 'bitops', 'chron', 'coda', 'colorspace', 'eRm',
'fExoticOptions', 'fMultivar', 'fOptions', 'Hmisc', 'ipred', 'lavaan', 'maptools',
'mathgraph', 'mgcv', 'mondate', 'MPV', 'mvtnorm', 'polspline', 'RANN',
'RcmdrPlugin.qual', 'rcom', 'RcppArmadillo', 'rgdal', 'rgeos', 'rgl',
'RGtk2Extras', 'rms', 'rpart', 'rscproxy', 'seriation', 'sets', 'sfsmisc', 'sp',
'spc', 'spdep', 'sudoku', 'tripack', 'TSA', 'tseries', 'tweedie', 'WriteXLS',
'XLConnect', 'XML', 'xts', 'Zelig'
Update all/some/none? [a/s/n]:
a
Warning: package 'rscproxy' is in use and will not be installed
also installing the dependency 'lattice'

trying URL 'http://stat.ethz.ch/CRAN/bin/windows/contrib/2.14/lattice_0.20-10.zip'
Content type 'application/zip' length 710282 bytes (693 Kb)

Installing from source – if all else fails

(this means hand-installing other packages, avoid if possible!)

Bioconductor - RMassBank

www.bioconductor.org/packages/release/bioc/html/RMassBank.html

Details

biocViews	Bioinformatics , MassSpectrometry , Metabolomics , Software
Version	1.0.0
In Bioconductor since	BioC 2.11 (R-2.15)
License	Artistic-2.0
Depends	rcdk,yaml,mzR,methods
Imports	XML,RCurl
Suggests	gplots,RMassBankData,xcms
System Requirements	OpenBabel
URL	
Depends On Me	
Imports Me	
Suggests Me	RMassBankData

Package Downloads

Package Source	RMassBank 1.0.0.tar.gz
Windows Binary	RMassBank 1.0.0.zip (32- & 64-bit)
MacOS 10.5 (Leopard) binary	RMassBank 1.0.0.tgz
Package Downloads Report	Download Stats

Installing from source – if all else fails

(this means hand-installing other packages, avoid if possible!)

```
#If this doesn't work, try installing from source
# http://www.bioconductor.org/packages/release/bioc/html/RMassBank.html
# download Package : > install.packages("C:/DATA/RMassBank/RMassBank_Versions/RMassBank_1.0.0.tar.gz",
# http://www.bioconductor.org/packages/release/bioc/html/RMassBank.html repos=NULL, type="source")
install.packages("RMassBank") Installing package(s) into 'C:/Users/schymaem/Documents/R/win-library/2.14'
                                (as 'lib' is unspecified)
                                * installing *source* package 'RMassBank' ...
                                ** R
                                ** inst
                                ** preparing package for lazy loading
                                Warning: package 'rcdk' was built under R version 2.14.2
                                Warning: package 'rcdklibs' was built under R version 2.14.2
                                Warning: package 'fingerprint' was built under R version 2.14.2
                                Warning: package 'iterators' was built under R version 2.14.2
                                Warning: package 'yaml' was built under R version 2.14.2
                                Warning: package 'Rcpp' was built under R version 2.14.2
                                ** help
                                *** installing help indices
                                ** building package indices ...
                                *** tangling vignette sources ...
                                'RMassBank.Rnw'
                                ** testing if installed package can be loaded
                                Warning messages:
                                1: package 'rcdk' was built under R version 2.14.2
                                2: package 'rcdklibs' was built under R version 2.14.2
                                3: package 'fingerprint' was built under R version 2.14.2
                                4: package 'iterators' was built under R version 2.14.2
                                5: package 'yaml' was built under R version 2.14.2
                                6: package 'Rcpp' was built under R version 2.14.2

                                * DONE (RMassBank)
                                > |
```

Loading RMassBank

```
30 # Once you have got through all that, you need to load the packages:
31 library("RMassBank")
32 library("RMassBankData")
33
```

32:25 (Top Level) ↕

R Script ↕

Console ~/ ↕

```
> # Once you have got through all that, you need to load the packages:
```

```
> library("RMassBank")
```

```
Loading required package: mzR
Loading required package: Rcpp
Loading required package: XML
Loading required package: RCurl
Loading required package: bitops
Loading required package: zoo
```

```
Attaching package: 'zoo'
```

```
The following object(s) are masked from 'package:base':
```

```
as.Date, as.Date.numeric
```

```
Loading required package:
Loading required package:
Loading required package:
Loading required package:
Loading required package:
```

```
Loading required package: iterators
```

```
Warning messages:
```

```
1: package 'Rcpp' was built under R version 2.14.2
2: package 'XML' was built under R version 2.14.2
3: package 'bitops' was built under R version 2.14.2
4: package 'zoo' was built under R version 2.14.2
5: package 'rcdk' was built under R version 2.14.2
6: package 'rcdklibs' was built under R version 2.14.2
7: package 'fingerprint' was built under R version 2.14.2
8: package 'iterators' was built under R version 2.14.2
```

```
> library("RMassBankData")
```


```
Warning message:
```

```
package 'RMassBankData' was built under R version 2.14.2
```

```
>
```


Finding out more about RMassBank

vignette(“RMassBank”)



RMassBank: The workflow by example

Michael Stravs

October 2, 2012

We'll kind-of work through this today...

1 Introduction	2
2 Installation and loading	2
3 Input files	3
3.1 LC/MS data	3
3.2 Compound list	4
3.3 Settings	4
4 The workflow	7
4.1 Mass spectrometry workflow	7
4.2 MassBank record workflow	14
5 Session information	16

Before getting started:

Install / Check OpenBabel

If you think you already have this, check e.g. (for Windows....)

- All Programs => OpenBabel 2.X.X
- C:\Program Files\OpenBabel-2.3.0

If you don't have this, please download (not *required* but it is *recommended*)

- <http://openbabel.org/>


The screenshot shows the OpenBabel website in a browser. The address bar displays 'openbabel.org/wiki/Main_Page'. The page title is 'Open Babel: The Open Source Chemistry Toolbox'. A navigation menu on the left includes links like 'Main Page', 'Get Open Babel', and 'Need Help?'. The main content area features a list of bullet points describing the software's capabilities and usage statistics. A large green arrow pointing downwards, labeled 'Download', is highlighted with a red rectangular border. To the right of the arrow are links for 'Official User Guide', 'Online docs', 'Buy Book', and 'Read Paper'. Further right is a diagram showing 'OBBase' connected to 'OBMol' and 'OE', with a 'Browse the API' link below it.

Before getting started:

Install / Check OpenBabel



Why OpenBabel?

- Convert SMILES (c1ccccc1) to structure
- C:\Program Files\OpenBabel-2.3.0



NOTE:

- By default, RMassBank will work without OpenBabel and has a back-up solution using CACTUS
- But, structures are less visually appealing



Before getting started:

Install / Check Proteowizard













Why Proteowizard?

- Convert raw data files to an open format (we chose mzML)
- Is one of most established and supported converters

**NOTE: If you will only ever use RMassBank once (i.e. today)
Don't bother installing! The test data is already converted!**

Check to see if you have it

- E.g. for me: C:\DATA\Program_Info\pwiz
- You want "MSConvertGUI.exe"

Name	Date modified	Type	Size
 MassLynxRaw.dll	27.01.2012 03:00	Application extens...	125 KB
 MassSpecDataReader.dll	27.01.2012 03:00	Application extens...	44 KB
 MathNet.Iridium.dll	26.04.2011 23:15	Application extens...	240 KB
 msaccess.exe	27.01.2012 03:01	Application	5'792 KB
 msbenchmark.exe	27.01.2012 03:03	Application	4'044 KB
 mscat.exe	27.01.2012 03:02	Application	4'025 KB
 msconvert.exe	27.01.2012 03:01	Application	5'755 KB
 MSConvertGUI.exe	27.01.2012 03:01	Application	54 KB
 MSConvertGUI.exe.manifest	09.01.2012 22:44	MANIFEST File	3 KB
 msdiff.exe	27.01.2012 03:01	Application	4'206 KB
 msdir.exe	27.01.2012 03:01	Application	4'166 KB
 MSFileReader.XRawfile2.dll	13.12.2011 06:22	Application extens...	237 KB

Before getting started:

Install / Check Proteowizard

No Proteowizard (pwiz)?: <http://proteowizard.sourceforge.net/downloads.shtml>

proteowizard.sourceforge.net/downloads.shtml

info download user docs dev docs contact

ProteoWizard Downloads

[Get ProteoWizard \(3 steps\)](#)

Step 1, Choose Download Type:

Windows installer (includes vendor reader support)

Step 2, Read License Agreements

Pwiz Core AB Sciex WiffReader Agilent MHDAC Bruker CompassXtract Thermo msFileReader Waters WRDAC

ADDENDUM TO APACHE LICENSE

To the best of our ability we deliver this software to you under the Apache 2.0 License listed below (the source code is available in the ProteoWizard project). This software does, however, depend on other software libraries which place further restrictions on its use and redistribution. By accepting the license terms for this software, you agree to comply with the restrictions imposed on you by the license agreements of the software libraries on which it depends:

Agilent Mass Hunter Data Access Component Library
Waters Raw Data Access Component Library
Bruker CompassXtract
Thermo-Scientific MSFileReader Library
Waters Raw Data Access Component Library
AB Sciex WIFF Reader Library

NOTE: If you do not plan to redistribute this software yourself, then you are the "end-user" in the above agreements.

Apache License - Covers Core ProteoWizard Tools and Library


Step 3, Click below to agree to Licenses and get ProteoWizard

I agree to the licensing terms, download ProteoWizard

Before getting started:

Install / Check Proteowizard


To use profile data (RMassBank deprofiler):



Before getting started:

Install / Check Proteowizard

To use centroid data (Pwiz deprofiler – example data prepared this way):



Before getting started

File names...

The files are used to identify a compound in the workflow

- Currently, we require a format like:

xxxxxxxxxxxxx_1234_.mzML

- _1234_ is the compound ID which is also in the compound list
- mzML is the format currently read by RMassBank
- Advanced / alternative formats possible
 - Provided in the GUI; more details from Steffen

Before getting started:

The compound list

Go back to the R / RStudio window

- Choose a working directory (change this!!! and create a new one if you like)
- Copy the compound list into this directory:

```
39 # Save a copy of the compound list|
40 # First, set a working directory (I make sure this exists)
41 setwd("C:/DATA/RMassBank/RMassBank_Demo")
42 # then copy the file there
43 file.copy(system.file("list/NarcoticsDataset.csv",
44                       package="RMassBankData"), "./Compoundlist.csv")
45
```

```
> setwd("C:/DATA/RMassBank/RMassBank_Demo")
> # then copy the file there
> file.copy(system.file("list/NarcoticsDataset.csv",
+                       package="RMassBankData"), "./Compoundlist.csv")
[1] TRUE
>
```


Before getting started:


The settings file (lucky last bit!!)

A number of settings are necessary for successful generation:

- Data processing
- Correct annotation of spectra

```
# Now, generate a settings file:  
RmbSettingsTemplate("mysettings.ini")
```

- Open up in any text editor and edit...



```
mysettings.ini - Notepad  
File Edit Format View Help  
# Sample configuration file for RMassBank.  
# Adapt this file to your needs.  
# NOTE: Do not indent with TAB characters! Use only spaces.  
# (If your editor converts TAB to a certain number of spaces, it's OK.)  
# Use a space after the colon.  
# Deprofile input data?  
# Leave empty if input data is already in "centroid" mode.  
# Use values deprofile.spline, deprofile.fwhm or deprofile.localMax to convert the  
input data with the  
# corresponding algorithm. See ?deprofile  
deprofile:  
# Deviation (in minutes) allowed the for retention time  
rtMargin: 0.4
```

Before getting started:

The settings file (lucky last bit!!)

```

mysettings.ini - Notepad
File Edit Format View Help

# Deviation (in minutes) allowed the for retention time
rtMargin: 0.4
# Systematic retention time shift
rtShift: -0.3

# Directory to openBabel. Required for creating molfiles for MassBank export.
# If no OpenBabel directory is given, RMassBank will attempt to use the CACTUS
webservice
# for SDF generation. You really should install openBabel though; the CACTUS
structures
# have explicit hydrogen atoms...
# Points to the directory where babel.exe (or the Linux "babel" binary) is located.
babeldir: |
# Example:
# babeldir: "C:\Program Files (x86)\OpenBabel-2.3.1\"

# which MassBank record version to use; version 2 is advised.
use_version: 2

# Include reanalyzed peaks?
use_rean_peaks: TRUE

# annotate the spectra files with (putative) molecular formulas for fragments?
add_annotation: TRUE

# Annotations for the spectrum:
annotations:
# Author etc. annotation
authors: Nomen Nescio, The Unseen University
copyright: Copyright (c) xxx
license: CC BY-SA
instrument: LTQ Orbitrap XL Thermo Scientific
instrument_type: LC-ESI-ITFT
confidence_comment: standard compound

```

These settings are OK for test data

Add the OpenBabel location here
CAREFUL – ‘ and “ as in example

Define record format (keep 2)

These are either TRUE or FALSE

Reanalyzed = formula + N₂O

Annotation adds formula to records

Change author and copyright to
something relevant for you...

Before getting started:

The settings file (lucky last bit!!)

```

mysettings.ini - Notepad
File Edit Format View Help
annotations:
  # Author etc. annotation
  authors: Nomen Nescio, The Unseen University
  copyright: Copyright (c) XXX
  license: CC BY-SA
  instrument: LTQ Orbitrap XL Thermo Scientific
  instrument_type: LC-ESI-ITFT
  confidence_comment: standard compound
  compound_class: N/A; Environmental Standard
  internal_id_fieldname: INTERNAL_ID
#
# HPLC annotations:
#
# example: lc_gradient: 90/10 at 0 min, 50/50 at 4 min
# example: lc_flow: 200 uL/min
# example: lc_solvent_a: water with 0.1% formic acid
# example: lc_column: XBridge C18 3.5um, 2.1x50mm, Waters
lc_gradient:
lc_flow:
lc_solvent_a:
lc_solvent_b:
lc_column:
# Prefix for MassBank accession IDs
entry_prefix: XX
ms_type: MS2
ionization: ESI
ms_dataprocessing:
  RECALIBRATE: loess on assigned fragments and MS1
# List of data-dependent scans in their order (relative to the parent scan), for
# annotation of the MassBank records
# For every data-dependent scan event, specify an element with:
# mode: fragmentation mode, e.g. CID
# ces: "short" format collision energy (for record title)

```

License: change for own records?
These settings are OK for test data

These are currently empty
You can use the examples for the example data. Solvent B was methanol!
These values will need adjusting for own chromatographic program

Choose your record letters!
E.g. Eawag = EA

These are OK for test data

Before getting started:

The settings file (lucky last bit!!)

```
mysettings.ini - Notepad
File Edit Format View Help
# List of data-dependent scans in their order (relative to the parent scan), for
# annotation of the MassBank records
# For every data-dependent scan event, specify an element with:
# mode: fragmentation mode, e.g. CID
# ces: "short" format collision energy (for record title)
# ce: "long" format collision energy (for annotation field)
# res: FT resolution
spectralList:
# First scan: CID 35% NCE, resolution 7500
- mode: CID
  ces: 35%
  ce: 35 % (nominal)
  res: 7500
# Second scan: HCD 15% NCE, resolution 7500
- mode: HCD
  ces: 15%
  ce: 15 % (nominal)
  res: 7500
# Third scan, etc.
- mode: HCD
  ces: 30%
  ce: 30 % (nominal)
  res: 7500
- mode: HCD
  ces: 45%
  ce: 45 % (nominal)
  res: 7500
- mode: HCD
  ces: 60%
  ce: 60 % (nominal)
  res: 7500
- mode: HCD
  ces: 75%
  ce: 75 % (nominal)
```

These settings are OK for test data

However:
This will need adjustment for own
experiments!

Before getting started:

The settings file (lucky last bit!!)

```

mysettings.ini - Notepad
File Edit Format View Help
# Shifts of the starting points for RMassBank accession numbers.
# Change these if you measure different adducts
accessionNumbersShifts:
  pH: 0 # [M+H]+: Accession numbers 1-14
  pM: 16 # [M]+: 17-30
  pNa: 32 # [M+Na]+: 33-46
  mH: 50 # [M-H]-: 51-64
  mFA: 66 # [M+FA]-: 67-80

# A list of known electronic noise peaks
electronicNoise:
- 189.825
- 201.725
- 196.875
# Exclusion width of electronic noise peaks (from unmatched peaks, or slightly varying)
# reanalysis)
electronicNoiseWidth: 0.3

# recalibration settings:
# recalibrate by: dppm or dmz
recalibrateBy: dppm

# recalibrate MS1:
# separately (separate)
# with common curve (common)
# do not recalibrate (none)
recalibrateMS1: common

# Custom recalibration function: You can overwrite the recal function by
# making any function which takes rcdata$recalfield ~ rcdata$mzFound.
# The settings define which recal function is used.
# Note: if recalibrateMS1 is "common", the setting "recalibrator: MS1" is
meaningless
# because the MS1 points will be recalibrated together with the MS2 points with
# the MS2 recalibration function.
recalibrator:
  MS1: recalibrate.loess
  MS2: recalibrate.loess

```

Fixed accession number: **XXYYYYZZ**;
XX = letters, **YYYY** = compound ID
ZZ = "shifts" – for us, these adducts

Remove consistent electronic noise
 Add the middle mass of the noise (often broad peaks, or slightly varying)
 Then add exclusion width –not too generous!


Recalibrate settings:
 Mass accuracy unit
 MS & MS/MS

These are OK for test data
 Changes necessary if recalibration needs adjusting

Before getting started:

The settings file (lucky last bit!!)

Finally: Save this before you go on!!!



```
# Define the multiplicity filtering level
# Default is 2 (peak occurs at least twice)
# Set this to 1 if you want to turn this option off.
# Set this to anything > 2 if you want harder filtering
multiplicityFilter: 2
```

Multiplicity Filtering settings e.g. require Peak to occur twice or more.

```
# Define the title format.
# You can use all entries from MassBank records as tokens
# plus the additional token RECORD_TITLE_CE, which is a shortened
# version of the collision energy specifically for use in the title.
# Every line is one entry and must have one token in curly brackets
# e.g. {CH$NAME} or {AC$MASS_SPECTROMETRY: MS_TYPE} plus optionally
# additional text in front or behind e.g.
# R={AC$MASS_SPECTROMETRY: RESOLUTION}
# If this is not specified, it defaults to a title of the format
# "Dinotefuran; LC-ESI-QFT; MS2; CE: 35%; R=35000; [M+H]+"
# Note how everything must be in "" here because otherwise the : are get
titleFormat:
- "{CH$NAME}"
- "{AC$INSTRUMENT_TYPE}"
- "{AC$MASS_SPECTROMETRY: MS_TYPE}"
- "CE: {RECORD_TITLE_CE}"
- "R={AC$MASS_SPECTROMETRY: RESOLUTION}"
- "{MS$FOCUSED_ION: PRECURSOR_TYPE}"
```

Change the title format
If you wish to display different
information in the title

```
# Define filter settings.
# For Orbitrap, settings of 15 ppm in low mass range, 10 ppm in high
# mass range, m/z = 120 as mass range division and 5 ppm for recalibrated
# data overall are recommended.
filtersettings:
  ppmHighMass: 10
  ppmLowMass: 15
  massRangeDivision: 120
  ppmFine: 5
  prelimCut: 1e4
  prelimCutRatio: 0
  fineCut: 0
  fineCutRatio: 0
  specokLimit: 1e4
  dbeminLimit: -0.5
```

These are all advanced settings
and do not need to be changed

Before getting started:

The settings file (lucky last bit!!)

Finally: Save this before you go on!!!

```

mysettings_ES.ini - Notepad
File Edit Format View Help

# Note how everything must be in "" here because otherwise the : are getting mangled!
titleFormat:
- "{CH$NAME}"
- "{AC$INSTRUMENT_TYPE}"
- "{AC$MASS_SPECTROMETRY: MS_TYPE}"
- "CE: {RECORD_TITLE_CE}"
- "R={AC$MASS_SPECTROMETRY: RESOLUTION}"
- "{MS$FOCUSED_ION: PRECURSOR_TYPE}"

# Define filter settings.
# For Orbitrap, settings of 15 ppm in low mass range, 10 ppm in high
# mass range, m/z = 120 as mass range division and 5 ppm for recalibrated
# data overall are recommended.
filtersettings:
  ppmHighMass: 10
  ppmLowMass: 15
  massRangeDivision: 120
  ppmFine: 5
  prelimCut: 1e4
  prelimCutRatio: 0
  fineCut: 0
  fineCutRatio: 0
  specokLimit: 1e4
  dbeminLimit: -0.5
  satelliteMzLimit: 0.5
  satelliteIntLimit: 0.05


# Define raw MS retrieval settings.
findMSRawSettings:
  ppmFine: 10
  mzCoarse: 0.5
  # fillPrecursorScan is FALSE for "good" mzML files which have all the info needed.
  # However, for example AB sciex files will have missing precursor scan information,
  # in which case fillPrecursorScan = TRUE is needed. Try it out.
  fillPrecursorScan: FALSE
  
```

These are all advanced settings and do not need to be changed

Adjust the boundaries for recalibration

**Change the MS retrieval
May be needed for different
Instruments than the default.**

RMassBank Workflow – Simple Form




Mass Spectrometry Workflow

Now we can get started!

Go back to your R session and “RMassBank_Demo.R”

- Select and run this whole snippet – but change the file name first!



```
52 # Start the MS Workflow
53 loadRmbSettings("mysettings_ES.ini")
54 # Start a new workspace
55 w <- newMsmsWorkspace()
56 # and load the files from the test data set
57 files <- list.files("test_data", pattern = "*.mzML")
58
59 # if you want to use the test data set
60 basename(files)
61 # and add the files to the workspace
62 w@files <- files
63 # if you want to use the test data set
64 #w@files <- files
65 # Then load the files
66 loadList(w@files)
67
```

```
Console C:/DATA/RMassBank/RMassBank_Demo/
> basename(files)
[1] "1_3_Chlorophenyl_piperazin_2818_pos.mzML"
[2] "1_3-Trifluoromethylphenyl_piperazin_2819_pos.mzML"
[3] "1_Benzylpiperazin_2820_pos.mzML"
[4] "Amitriptylin_2821_pos.mzML"
[5] "Amphetamin_2822_pos.mzML"
[6] "Benzoylecgonin_2823_pos.mzML"
[7] "Cocain_2817_pos.mzML"
[8] "Dextromethorphan_2824_pos.mzML"
[9] "EDDP_2_Ethyl_1_5_dimethyl_3_3_diphenylpyrrolinium_2825_pos.mzML"
[10] "Ephedrin_2758_pos.mzML"
[11] "Ketamin_2826_pos.mzML"
[12] "Mephedron_4_Methylmethcathinon_2827_pos.mzML"
[13] "Methadon_2828_pos.mzML"
[14] "Methamphetamin_2829_pos.mzML"
[15] "Naltrexon_2830_pos.mzML"
> # and add the files to the workspace
> w@files <- files
```

Mass Spectrometry Workflow

msmsWorkflow: 8 Steps in summary

- Step 1: Search all `files` for MS/MS spectra of compound
- Step 2: Molecular formula fitting for all peaks (large tolerance)
- Step 3: Analyzed spectra aggregated into a list
- Step 4: Recalibration curve is calculated; all spectra are recalibrated
- Step 5: Recalibrated spectra are reanalyzed (new formula fit post calibration)
- Step 6: Reanalyzed spectra aggregated. Electronic noise removed.
 - Back-up saved as `archiveName.RData`
- Step 7: Fail peaks are reanalyzed (formula + 2N, O)
 - Back-up saved as `archiveName_RA.RData`
- Step 8: Multiplicity Filtering: peaks occurring once only are removed.
 - Back-up saved as `archiveName_RF.RData`
 - High Intensity Fail Peaks: `archiveName_Failpeaks.csv`

Mass Spectrometry Workflow

8 steps in RMassBank code...

- Steps 1:4 – up to recalibration

```
# Start the workflow with steps 1 to 4:
w <- msmsWorkflow(w, mode="pH", steps=c(1:4), archivename = "pH_narcotics")
```

- Warning: this can take a while....

- Progress:


```
> # Start the workflow with steps 1 to 4:
> w <- msmsWorkflow(w, mode="pH", steps=c(1:4), archivename = "pH_narcotics")
msmsWorkflow: Step 1
|
| 0%3: >C:/Users/schymaem/Documents/R/win-
library/2.14/RMassBankData/spectra/1_3_Chlorophenyl_piperazin_2818_pos.mzML
|
| 7%3: >C:/Users/schymaem/Documents/R/win-
library/2.14/RMassBankData/spectra/1_3-Trifluoromethylphenyl_piperazin_2819_pos.mzML
|
| 13%3: >C:/Users/schymaem/Documents/R/win-
library/2.14/RMassBankData/spectra/1_Benzylpiperazin_2820_pos.mzML
|
|=====| 20%3:
>C:/Users/schymaem/Documents/R/win-library/2.14/RMassBankData/spectra/Amitriptylin_2821
```

Mass Spectrometry Workflow


Steps 1:4 of RMassBank code...

- Recalibration plot at the end


MS2 scatterplot




MS2 density



MS1 scatterplot



MS1 density



Mass Spectrometry Workflow

Steps 5:8 in RMassBank code...

- Steps 5:8 – complete the MS workflow

```
# Keep going with steps 5-8:
w <- msmsWorkflow(w, mode="pH", steps=c(5:8), archivename = "pH_narcotics")
```

- Warning: this can also take a while....

- Progress:


```
> # Keep going with steps 5-8:
> w <- msmsWorkflow(w, mode="pH", steps=c(5:8), archivename = "pH_narcotics")
msmsWorkflow: Step 5
|
| 0% |
| 7% |
|====|
|====|
|====|
|====|
|====|
|====|
|====|
|====|
|====|
|====|
```

13%
20%
27%
33%
40%
47%
53%
60%

- Now we need to look at the “Fail Peaks”

Mass Spectrometry Workflow

File listing should now look something like:



- o Open up the fail peaks (outside R) and have a look

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	
1	OK	name	cpdID	scan	mzFound	formula	mzCalc	dppm	dbe	mz	int	formulaC	parentSca	aMax	mzCenter	
2		407	2817	558	182.1191	NA	NA	NA	NA	182.1193	117786.4	0	556	9902212	304.1543	
3		2231	2821	963	233.1351	NA	NA	NA	NA	233.1355	255822.6	0	961	16383558	278.1903	
4		3491	2824	723	215.1468	NA	NA	NA	NA	215.1471	242559.4	0	721	15875343	272.2009	
5		3988	2825	725	249.1548	NA	NA	NA	NA	249.1553	516037.1	0	723	23457888	278.1903	
6		4110	2825	730	193.0993	NA	NA	NA	NA	193.0995	138038.3	0	723	5673618	278.1903	
7		5543	2828	938	265.1612	NA	NA	NA	NA	265.1618	146493.5	0	933	4158767	310.2165	
8		6503	2830	352	211.0738	C12H9N3O	211.074	-0.9633		13	211.0741	24055.95	1	346	195149.4	342.17

Mass Spectrometry Workflow

Fail Peaks



Compound ID MS/MS Scan m/z & I of MS/MS peak I, m/z of precursor

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1	OK	name	cpdID	scan	mzFound	formula	mzCalc	dppm	dbe	mz	int	formulaC	parentScan	aMax	mzCenter
2		407	2817	558	182.1191	NA	NA	NA	NA	182.1193	117786.4	0	556	9902212	304.1543
3		2231	2821	963	233.1351	NA	NA	NA	NA	233.1355	255822.6	0	961	16383558	278.1903
4		3491	2824	723	215.1468	NA	NA	NA	NA	215.1471	242559.4	0	721	15875343	272.2009
5		3988	2825	725	249.1544	NA	NA	NA	NA	249.1553	516837.1	0	723	23457888	278.1903
6		4110	2825	730	193.0993	NA	NA	NA	NA	193.0995	138038.3	0	723	5673618	278.1903
7		5543	2828	938	265.1617	NA	NA	NA	NA	265.1618	146493.5	0	933	4158767	310.2165
8		6503	2830	352	211.0738	C12H9N3O	211.074	-0.9633	13	211.0741	24055.95	1	346	195149.4	342.17

No formula calculated; not replicated (only present once) and low I compared with parent => "true" fail peaks
Look at raw data of "2825"

Mass Spectrometry Workflow


Fail Peak



Mass Spectrometry Workflow

Fail Peaks

EDDP_2_Ethyl_1_5_dimethyl_3_3_diphenylpyrrolinium_2825_pos #730 RT: 6.59 AV: 1 NL: 5.59E6
T: FTMS + p ESI d Full ms2 278.19@hcd90.00 [50.00-290.00]



“Fail peak” is very minor part of spectrum, already within noise levels

Mass Spectrometry Workflow


Fail Peaks

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1	OK	name	cpdID	scan	mzFound	formula	mzCalc	dppm	dbe	mz	int	formulaC	parentSca	aMax	mzCenter
2		407	2817	558	182.1191	NA	NA	NA	NA	182.1193	117786.4	0	556	9902212	304.1543
3		2231	2821	963	233.1351	NA	NA	NA	NA	233.1355	255822.6	0	961	16383558	278.1903
4		3491	2824	723	215.1468	NA	NA	NA	NA	215.1471	242559.4	0	721	15875343	272.2009
5		3988	2825	725	249.1548	NA	NA	NA	NA	249.1553	516037.1	0	723	23457888	278.1903
6		4110	2825	730	193.0993	NA	NA	NA	NA	193.0995	138038.3	0	723	5673618	278.1903
7		5543	2828	938	205.1012	NA	NA	NA	NA	205.1018	140493.5	0	933	4158767	310.2165
8		6503	2830	352	211.0738	C12H9N3O	211.074	-0.9633	13	211.0741	24055.95	1	346	195149.4	342.17


Formula calculated, but also low intensity compared with precursor
Go to raw data...

Mass Spectrometry Workflow

Fail Peak
(with formula)




Naltrexon_2830_pos #352 RT: 3.19 AV: 1 NL: 1.95E5
T: FTMS + p ESI d Full ms2 342.17@hcd75.00 [50.00-355.00]




Mass Spectrometry Workflow

Fail peak (with formula)

Naltrexon_2830_pos #352 RT: 3.19 AV: 1 NL: 1.95E5
T: FTMS + p ESI d Full ms2 342.17@hcd75.00 [50.00-355.00]



Naltrexon_2830_pos #353 RT: 3.19 AV: 1 NL: 2.06E5
T: FTMS + p ESI d Full ms2 342.17@hcd90.00 [50.00-355.00]



Mass Spectrometry Workflow

Fail Peaks


Include this “fail peak” by entering a 1 in the OK column

- Save under a new name...
- e.g. “pH_narcotics_Failpeaks_wOKs.csv”

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O
1	OK	name	cpdID	scan	mzFound	formula	mzCalc	dppm	dbe	mz	int	formulaC	parentSca	aMax	mzCe
2	0	407	2817	558	182.1191	NA	NA	NA	NA	182.1193	117786.4	0	556	9902212	304.
3	0	2231	2821	963	233.1351	NA	NA	NA	NA	233.1355	255822.6	0	961	16383558	278.
4	0	3491	2824	723	215.1468	NA	NA	NA	NA	215.1471	242559.4	0	721	15875343	272.
5	0	3988	2825	725	249.1548	NA	NA	NA	NA	249.1553	516037.1	0	723	23457888	278.
6	0	4110	2825	730	193.0993	NA	NA	NA	NA	193.0995	138038.3	0	723	5673618	278.
7	0	5543	2828	938	265.1612	NA	NA	NA	NA	265.1618	146493.5	0	933	4158767	310.
8	1	6503	2830	352	211.0738	C12H9N3C	211.074	-0.9633	13	211.0741	24055.95	1	346	195149.4	34
9															

We are now ready to go onto the Record Creation part of the workflow...

RMassBank Workflow – Simple Form



MassBank Record Workflow

Start a new MB workspace and reset the “Infolists” (compound information)

- Reset and load “Infolists” (compound information) – we will use a part-filled one

```
79
80 # Next, start the MassBank Record Workflow:
81 mb <- newMbWorkspace(w)
82 mb <- resetInfolists(mb)
83 # To speed up this example, we have provided a partially complete list
84 mb <- loadInfolists(mb, system.file("infolists_incomplete",
85                                     package="RMassBankData"))
```

- Add fail peaks, if applicable, and start the workflow:

```
90 # now add the fail peaks from above, if you have included any
91 mb <- addPeaks(mb, "pH_narcotics_Failpeaks_wOKs.csv")
92
93 # Start the record generation workflow
94 mb <- mbWorkflow(mb, infolist_path="./Narcotics_infolist.csv")
95 |
```

MassBank Record Workflow

Progress:

```

Console C:/DATA/RMassBank/RMassBank_Demo/
3
3
2824: smiles
3
3
2825: dbname
3
3
3
3
2758: smiles
3
3
2826: smiles
mbWorkflow: Step 2
The file ./Narcotics_infolist.csv was generated with new compound information. Please
check and edit the table, and add it to your infolist folder.

```

This shows the source of the data for querying CACTUS and CTS: Most come from SMILES, but the back-up of a “name” (dbname) is useful in some cases

There is new information to check manually
Please open the table in Excel/OpenOffice...

MassBank Record Workflow

Infolist Editing (1) - Names

	A	B	C	D	E	F	G	H	I	J	K	L	M
1		id	dbcas	dbname	dataused	COMMENT	COMMENT	CH\$NAME	CH\$NAME	CH\$NAME	CH\$COMP	CH\$FORM	CH\$EXACT
2	1	2819	15532-75-1	1-(3-Triflu	smiles	standard c	2819	1-[3-(triflu	1-(3-(Trifluoromethyl	N/A; Envir	C11H13F3	230.1031	
3	2	2821	50-48-6	Amitriptyl	smiles	standard c	2821			N/A; Envir	C20H23N	277.183	
4	3	2822	300-62-9	Amphetar	smiles	standard c	2822	Amphetar	Amfetami (1-methyl	N/A; Envir	C9H13N	135.1048	
5	4	2823	519-09-5	Benzoylec	dbname	standard c	2823	Benzoylec	None	3-benzoy	N/A; Envir	C16H19NC	289.1314
6	5	2817	50-36-2	Cocaine	dbname	standard c	2817	Cocaine	Unknown	(1S,3S,4R,	N/A; Envir	C17H21NC	303.1471
7	6	2824	125-71-3	Dextrome	smiles	standard c	2824			N/A; Envir	C18H25NC	271.1936	
8	7	2825	30223-73-	EDDP (2-E	dbname	standard c	2825			N/A; Envir	C20H23N	277.183	
9	8	2758	299-42-3	Ephedrine	smiles	standard c	2758	Ephedrine	(1R,2S)-2-(methylam	N/A; Envir	C10H15NC	165.1154	
10	9	2826	6740-88-1	Ketamine	smiles	standard c	2826			N/A; Envir	C13H16Cl	237.092	
11													

Fix up the names here

- “None”, “Unknown” unnecessary
- At least one name should be present (copy from dbname)

H	I	J	K
CH\$NAME	CH\$NAME	CH\$NAME	CH\$COMP
1-[3-(triflu	1-(3-(Trifluoromethyl	N/A; Envir	
Amitriptyline		N/A; Envir	
Amphetar	Amfetami (1-methyl	N/A; Envir	
Benzoylec	3-benzoyloxy-8-me	N/A; Envir	
Cocaine	(1S,3S,4R,5R)-3-ben	N/A; Envir	
Dextromethorphan		N/A; Envir	
EDDP (2-Ethylidene-1,5-dimet		N/A; Envir	
Ephedrine	(1R,2S)-2-(methylam	N/A; Envir	
Ketamine		N/A; Envir	

MassBank Record Workflow

Infolist Editing (2) – double-check CAS Numbers (if they were in original csv)

C	D	E	F	G	H	I	J	K	L	M	N	O	P
dbcas	dbname	dataused	COMMENT	COMMENT	CH\$NAME	CH\$NAME	CH\$NAME	CH\$COMP	CH\$FORM	CH\$EXACT	CH\$SMILE	CH\$IUPAC	CH\$LINK.C
15532-75-9	1-(3-Trifluoromethylphenyl)ethanone	smiles	standard	2819	1-[3-(trifluoromethyl)phenyl]ethanone	1-(3-(Trifluoromethyl)phenyl)ethanone		N/A; Envir	C11H13F3O	230.1031	c1c(C(F)(F)F)c(C)cc1	InChI=1S/C11H13F3O	15532-75-9
50-48-6	Amitriptyline	smiles	standard	2821	Amitriptyline			N/A; Envir	C20H23N	277.183	C1(\c2c(C)ccn2)CC1	InChI=1S/C20H23N	50-48-6
300-62-9	Amphetamine	smiles	standard	2822	Amphetamine	Amfetami (1-methylamino)ethanone		N/A; Envir	C9H13N	135.1048	c1(ccccc1)CN	InChI=1S/C9H13N	300-62-9
519-09-5	Benzoylcholine	dbname	standard	2823	Benzoylcholine	3-benzoyloxy-8-methyl-8-azabicyclo[3.2.1]octane		N/A; Envir	C16H19NO	289.1314	CN1C2CCC(C2)OC(=O)c1ccccc1	InChI=1S/C16H19NO	519-09-5
50-36-2	Cocaine	dbname	standard	2817	Cocaine	(1S,3S,4R,5R)-3-benzoyloxy-N-methyl-1,2,3,4-tetrahydroisoquinoline		N/A; Envir	C17H21NO	303.1471	CN1C2CCC(C2)OC(=O)c1ccccc1	InChI=1S/C17H21NO	50-36-2
125-71-3	Dextromethorphan	smiles	standard	2824	Dextromethorphan			N/A; Envir	C18H25NO	271.1936	c12[C@H](C)CCN1C2	InChI=1S/C18H25NO	125-71-3
30223-73-1	EDDP (2-Ethylidene-1,5-dimethylpiperazine)	dbname	standard	2825	EDDP (2-Ethylidene-1,5-dimethylpiperazine)			N/A; Envir	C20H23N	277.183	None	None	
299-42-3	Ephedrine	smiles	standard	2758	Ephedrine	(1R,2S)-2-(methylamino)propane-1,1-diol		N/A; Envir	C10H15NO	165.1154	c1([C@H](O)C)N(C)C1	InChI=1S/C10H15NO	299-42-3
6740-88-1	Ketamine	smiles	standard	2826	Ketamine			N/A; Envir	C13H16ClN	237.092	c1([C@@H](O)C)N(C)C1	InChI=1S/C13H16ClN	33643-46-8

MassBank Record Workflow

Infolist Editing (3) – reality-checking entries & delete “rubbish”

M	N	O	P	Q	R	S	T	U	V	W	X	Y	Z
CH\$EXACT	CH\$SMILE	CH\$IUPAC	CH\$LINK.C	CH\$LINK.C	CH\$LINK.F	CH\$LINK.F	CH\$LINK.I	CH\$LINK.F	CH\$LINK.F	CH\$LINK.F	CH\$LINK.F	CH\$LINK.F	CH\$LINK.F
230.1031	c1c(C(F)(F	InChI=1S/	15532-75-9										
277.183	C1(\c2c(C	InChI=1S/	50-48-6	2666		D07448							
135.1048	c1(ccccc1)	InChI=1S/	300-62-9	CHEBI	Delete this								
289.1314	CN1C2CCC	InChI=1S/	519-09-5										
303.1471	CN1C2CCC	InChI=1S/	50-36-2	27958		C01416							
271.1936	c12[C@]3	InChI=1S/	125-71-3			D03742							
277.183	None	None		None	None	None	None	None	None	None	None	None	None
165.1154	c1([C@H]	InChI=1S/	299-42-3	15407		C01575							
237.092	c1([C@@]	InChI=1S/	33643-46-8										


Check this - ChemSpider & PubChem ID usually similar

CHEBI

CID:2160 KRMDCWI 2075
CID:3007 KWTSXDU 13852819

MassBank Workflow

Infolist Editing (4) – checking “Amphetamine” in ChemSpider



RSC Synt


RSC

[About](#) |
 [More Searches](#) |
 [Web APIs](#) |
 [Help](#)

Retrieved ChemSpider ID is correct, so no changes necessary

eg. Pyrid

Search term: **Amphetamine** (Found by approved synonym) ?



Amphetamine

ChemSpider ID: **13852819**
 Molecular Formula: C₉H₁₃N
 Average mass: 135.206207 Da
 Monoisotopic mass: 135.104797 Da

▼ Systematic name
 1-Phenyl-2-propanamine

- ▶ SMILES and InChIs
- ▶ Cite this record

? 2D 3D Save Zoom

0 of 1 defined stereocentres

Want to comment on this record?

Leave Feedback

MassBank Record Workflow

Infolist Editing (5) – fill in missing entries for EDDP - search by CAS

M	N	O	P	Q	R	S	T	U	V	W	X
CH\$EXACT	CH\$SMILE	CH\$IUPAC	CH\$LINK.C	CH\$LINK.C	CH\$LINK.F	CH\$LINK.F	CH\$LINK.L	CH\$LINK.F	CH\$LINK.I	CH\$LINK.CHEMSPIDER	
230.1031	c1c(C(F)(F	InChI=1S/15532-75-9						CID:4296	KKIMDKM	4145	
277.183	C1(\c2c(C	InChI=1S/150-48-6	2666		D07448			CID:2160	KRMDCWI	2075	
135.1048	c1(C(C(C(C	InChI=1S/300-63-8			D07445			CID:3007	KWTSXDU	13852819	
289.1314	CN1C2CC	InChI=1S/1519-09-5						CID:2340	GVGYEFKI	2250	
303.1471	CN1C2CC	InChI=1S/150-36-2	27958		C01416			CID:44622	ZPUCINDJ	10194104	
271.1936	c12[C@]3	InChI=1S/125-71-2			D03742			CID:69161	MKXZASY	13109865	
277.183	None	None	None	None	None	None	None	CID:None	None		
165.1154	c1([C@H]	InChI=1S/1299-42-3	15407		C01575			CID:9294	KWGRBVC	8935	
237.092	c1([C@@]	InChI=1S/133643-46-8							YQEZLKZA	158414	

Search ChemSpider and PubChem with CAS 30223-73-5 (in csv)

MassBank Record Workflow

Infolist Editing (5) – EDDP - search ChemSpider by CAS: 30223-73-5

www.chemsp

3 hits found.
Search terms: 30223-
(Found by synonym)


Grid Tile Table

ID Structure

4526936




- 0/1 defined



23254962




- 0/1 defined



23621123




- 0/1 defined



About More Searches Web APIs Help

eg. Pyridi



Pyrrolidine, 1,5-dimethyl-2-ethylidene-3,3-diphenyl

ChemSpider ID: 4526936

Molecular Formula: C₂₀H₂₃N
Average mass: 277.40329 Da
Monoisotopic mass: 277.1830

Systematic name
(2Z)-2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine

SMILES and InChIs

Cite this record

Want to comment on this record?

Leave Feedback

SMILES and InChIs

SMILES:

C/C=C\1/C(CC(N1C)C)(c2ccccc2)c3ccccc3

Copy

Std. InChI:

InChI=1S/C20H23N
/c1-4-19-20(15-16(2)21(19)3,17-11-7-5-8-12-1'
13-9-6-10-14-18/h4-14,16H,15H2,1-3H3

Copy

Std. InChIKey:

AJRJPORIQGYFMT-PVOVUMCXSA-N

Copy

Cite this record

Names and Identifiers

Names and Synonyms

Validated by Experts, Validated by Users, Non-Validated, Removed by Experts

EDDP


(2Z)-2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine

Pyrrolidine, 1,5-dimethyl-2-ethylidene-3,3-diphenyl

MassBank Record Workflow

=> <http://pubchem.ncbi.nlm.nih.gov/search/search.cgi#>

Infolist Editing (5) – EDDP - search PubChem by CAS: 30223-73-5



2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine - Compound



Summary (CID 5378015)

You can copy this straight from the URL!


Also known as: (Z)-2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, 30223-73-5, EDDP, Eddp-3,3, 2-Ethylidene-1,5-dimethyl-3,3-DPP, AC1NTC6G

Molecular Formula: C₂₀H₂₃N Molecular Weight: 277.40332 InChIKey: AJRJPORIQGYFMT-PVOVUMCXSA-N

Table of Contents [Show subcontent titles](#)

- [Identification](#)
- [Related Records](#)
- [Biomedical Effects and Toxicity](#)
- [Literature](#)
- [Classification](#)
- [Chemical and Physical Properties](#)

[Expand all sub-sections](#)



MassBank Record Workflow

Infolist Editing (5) – fill in missing entries for EDDP – now filled in?

135.1048	c1(ccccc1)	InChI=1S/	300-62-9	CHEBI		D07445		CID:3007	KWTSXDU	13852819
289.1314	CN1C2CCC	InChI=1S/	519-09-5					CID:2340	GVGYEFKI	2250
303.1471	CN1C2CCC	InChI=1S/	50-36-2	27958		C01416		CID:44622	ZPUCINDJ	10194104
271.1936	c12[C@]3	InChI=1S/	125-71-3			D03742		CID:69161	MKXZASY	13109865
277.183	None	None		None	None	None	None	CID:None	None	
165.1154	c1([C@H])	InChI=1S/	299-42-3	15407		C01575		CID:9294	KWGRBVC	8935
237.092	c1([C@@])	InChI=1S/	33643-46-8						YQEZLKZA	158414

	M	N	O	P	Q	R	S	T	U	V	W	X
	CH\$EXACT	CH\$SMILE	CH\$IUPAC	CH\$LINK.C	CH\$LINK.C	CH\$LINK.F	CH\$LINK.F	CH\$LINK.L	CH\$LINK.F	CH\$LINK.I	CH\$LINK.CHEMSPIDER	
	230.1031	c1c(C(F)(F	InChI=1S/	15532-75-9					CID:4296	KKIMDKM	4145	
	277.183	C1(\c2c(C	InChI=1S/	50-48-6	2666		D07448		CID:2160	KRMDCWI	2075	
	135.1048	c1(ccccc1)	InChI=1S/	300-62-9			D07445		CID:3007	KWTSXDU	13852819	
	289.1314	CN1C2CCC	InChI=1S/	519-09-5					CID:2340	GVGYEFKI	2250	
	303.1471	CN1C2CCC	InChI=1S/	50-36-2	27958		C01416		CID:44622	ZPUCINDJ	10194104	
	271.1936	c12[C@]3	InChI=1S/	125-71-3			D03742		CID:69161	MKXZASY	13109865	
	277.183	C/C=C\1/C	InChI=1S/	30223-73-5					CID:53780	AJRJPORIC	4526936	
	165.1154	c1([C@H])	InChI=1S/	299-42-3	15407		C01575		CID:9294	KWGRBVC	8935	
	237.092	c1([C@@])	InChI=1S/	33643-46-8						YQEZLKZA	158414	

MassBank Record Workflow

Infolist Editing (6) – Check missing PubChem entry?

M	N	O	P	Q	R	S	T	U	V	W	X
CH\$EXACT	CH\$SMILE	CH\$IUPAC	CH\$LINK.C	CH\$LINK.C	CH\$LINK.F	CH\$LINK.F	CH\$LINK.L	CH\$LINK.F	CH\$LINK.I	CH\$LINK.C	CHEMSPIDER
230.1031	c1c(C(F)(F	InChI=1S/	15532-75-9					CID:4296	KKIMDKM	4145	
277.183	C1(\c2c(C	InChI=1S/	50-48-6	2666		D07448		CID:2160	KRMDCWI	2075	
135.1048	c1(ccccc1)	InChI=1S/	300-62-9			D07445		CID:3007	KWTSXDU	13852819	
289.1314	CN1C2CCC	InChI=1S/	519-09-5					CID:2340	GVGYEFKI	2250	
303.1471	CN1C2CCC	InChI=1S/	50-36-2	27958		C01416		CID:44622	ZPUCINDJ	10194104	
271.1936	c12[C@]3	InChI=1S/	125-71-3			D03742		CID:69161	MKXZASY	13109865	
277.183	C/C=C\1/C	InChI=1S/	30223-73-5					CID:53780	AJRJPORIC	4526936	
165.1154	c1([C@H]	InChI=1S/	299-42-3	15407		C01575		CID:9294	KWGRBVC	8935	
237.092	c1([C@@]	InChI=1S/	33643-46-8						YQEZLKZA	158414	

pubchem.ncbi.nlm.nih.gov/summary/summary.cgi?cid=182137

NCBI

PubChem Compound

Search

Limits Advanced search

(S)-Ketamine - Compound Summary (CID 182137)

Also known as: Esketamine, (-)-Ketamine, (S)-(-)-Ketamine, l-Ketamine, CHEBI:60799, cyclohexanone, 2-(2-chlorophenyl)-2-(methylamino)-, (2S)-, 33643-46-8

Molecular Formula: C₁₃H₁₆ClNO Molecular Weight: 237.72524 InChIKey: YQEZLKZALYSWHR-ZDUSSCGKSA-N

Links and Relations

Properties
Compound ID: 182137

MassBank Record Workflow

Infolist editing finished! Save under a new name and restart workflow

- Moving Infolists into a separate Infolist folder is recommended!

```

100
101 mb <- resetInfolists(mb)
102 mb <- loadInfolists(mb, "C:/DATA/RMassBank/InfoLists_demo")
103 # NOTE: if you didn't want to correct the whole list, use this one:
104 #mb <- loadInfolists(mb, system.file("infolists", package="RMassBankData"))
105 # and start the workflow again - "no new data added" is a good sign!
106 mb <- mbWorkflow(mb)
107
> mb <- resetInfolists(mb)
> mb <- loadInfolists(mb, "C:/DATA/RMassBank/InfoLists_demo")
> # NOTE: if you didn't want to correct the whole list, use this one:
> #mb <- loadInfolists(mb, system.file("infolists", package="RMassBankData"))
> # and start the workflow again - "no new data added" is a good sign!
> mb <- mbWorkflow(mb)
mbWorkflow: Step 1
mbWorkflow: Step 2
No new data added.
mbWorkflow: Step 3
mbWorkflow: Step 4
Compiling: 1_3_Chlorophenyl_piperazin_2818_pos.mzML
Compiling: 1_3-Trifluoromethylphenyl_piperazin_2819_pos.mzML
Compiling: 1_Benzylpiperazin_2820_pos.mzML
Compiling: Amitriptylin_2821_pos.mzML
Compiling: Amphetamin_2822_pos.mzML
Compiling: Benzoylecgonin_2823_pos.mzML
3
3
3
mbWorkflow: Step 5
mbWorkflow: Step 6
mbWorkflow: Step 7
mbWorkflow: Step 8
>

```

And that's it – DONE!

MassBank Record Workflow

Take a look at your new records!

- You can open these in any text editor

The screenshot displays a Windows File Explorer window with the following structure:

- Root: C:\DATA\RMassBank\RMassBank_Demo
- Files in main window:
 - Compoundlist.csv (1.8 KB)
 - mysettings.ini (4.7 KB)
 - mysettings_ES.ini
 - pH_narcotics.RData
 - pH_narcotics_Failpeaks.csv
 - pH_narcotics_RA.RData
 - pH_narcotics_RF.RData
 - pH_narcotics_Failpeaks_wOK
 - Narcotics_infolist.csv (3.7 KB)
 - EA (Folder)
- Sub-window (EA folder):
 - moldata (Folder)
 - recdata (Folder)
- Detailed view of 'recdata' folder:

Name	Date modified	Type
EA275801.txt	20.11.2012 17:05	Text Document
EA275802.txt	20.11.2012 17:05	Text Document
EA275803.txt	20.11.2012 17:05	Text Document
EA275804.txt	20.11.2012 17:05	Text Document
EA275805.txt	20.11.2012 17:05	Text Document
EA275806.txt	20.11.2012 17:05	Text Document
EA275807.txt	20.11.2012 17:05	Text Document
EA275808.txt	20.11.2012 17:05	Text Document
EA275809.txt	20.11.2012 17:05	Text Document
EA275810.txt	20.11.2012 17:05	Text Document
EA275811.txt	20.11.2012 17:05	Text Document
EA275812.txt	20.11.2012 17:05	Text Document
EA275813.txt	20.11.2012 17:05	Text Document

MassBank Records “recdata”

```

EA275805.txt - Notepad
File Edit Format View Help
ACCESSION: EA275805
RECORD_TITLE: Ephedrine; LC-ESI-ITFT; MS2; 60%; R=7500; [M+H]+
DATE: 2012.11.20
AUTHORS: M. Stravs, E. Schymanski, H. Singer, Eawag
LICENSE: CC BY-SA
COPYRIGHT: Copyright (c) Eawag, 2012
COMMENT: CONFIDENCE standard compound
COMMENT: INTERNAL_ID 2758
CH$NAME: Ephedrine
CH$NAME: (1R,2S)-2-(methylamino)-1-phenyl-1-propanol
CH$COMPOUND_CLASS: N/A; Environmental standard
CH$FORMULA: C10H15NO
CH$EXACT_MASS: 165.1154
CH$SMILES: c1([C@H])([C@@H])(N)C(C)C1
CH$IUPAC: InChI=1S/C10H15NO/m0/s1
CH$LINK: CAS 299-42-3
CH$LINK: CHEBI 15407
CH$LINK: KEGG C01575
CH$LINK: PUBCHEM CID:9294
CH$LINK: INCHIKEY KWGRBVOPPL
CH$LINK: CHEMSPIDER 8935
AC$INSTRUMENT: LTQ Orbitrap
AC$INSTRUMENT_TYPE: LC-ESI-IT
AC$MASS_SPECTROMETRY: MS_TYP
AC$MASS_SPECTROMETRY: IONIZA
AC$MASS_SPECTROMETRY: ION_MO
AC$MASS_SPECTROMETRY: FRAGME
AC$MASS_SPECTROMETRY: COLLIS
AC$MASS_SPECTROMETRY: RESOLU
AC$CHROMATOGRAPHY: COLUMN_NA
AC$CHROMATOGRAPHY: FLOW_GRAD
5/95 at 25 min, 90/10 at 25.
AC$CHROMATOGRAPHY: FLOW_RATE
AC$CHROMATOGRAPHY: RETENTION
AC$CHROMATOGRAPHY: SOLVENT A
AC$CHROMATOGRAPHY: SOLVENT B
MS$FOCUSED_ION: BASE_PEAK 16
MS$FOCUSED_ION: PRECURSOR_M/Z 166.1226
MS$FOCUSED_ION: PRECURSOR_TYPE [M+H]+
MS$DATA_PROCESSING: RECALIBRATE loess on assigned fragments and MS1
MS$DATA_PROCESSING: WHOLE RmassBank
PK$ANNOTATION: m/z num {formula mass error(ppm)}
56.0495 1 C3H6N+ 56.0495 -0.46
70.0652 1 C4H8N+ 70.0651 0.49
91.0542 1 C7H7+ 91.0542 -0.51
93.07 1 C7H9+ 93.0699 1.22
115.0542 1 C9H7+ 115.0542 -0.58
117.0699 1 C9H9+ 117.0699 0.11
118.0647 1 C8H8N+ 118.0651 -3.86
132.081 1 C9H10N+ 132.0808 1.55
133.0886 1 C9H11N+ 133.0886 0.14
135.0804 1 C9H11O+ 135.0804 -0.68
148.1121 1 C10H14N+ 148.1121 -0.04
PK$NUM_PEAK: 11
PK$PEAK: m/z int. rel.int.
56.0495 151421.4 24
70.0652 152358.5 24
91.0542 150051.6 24
93.07 24399.9 3
115.0542 229068.4 37
117.0699 1548461.5 251
118.0647 12712.5 2
132.081 35161.8 5
133.0886 1423545.6 231
135.0804 94627.2 15
148.1121 6145796.5 999
//

```


MassBank Record Structures “moldata”

RMassBank_Demo ▶ EA ▶ moldata

Search moldata

Share with ▾ Burn New folder

Name	Date modified	Type
2758.mol	20.11.2012 17:05	MOL File
2817.mol	20.11.2012 17:05	MOL File
2818.mol	20.11.2012 17:05	MOL File
2819.mol	20.11.2012 17:05	MOL File
2820.mol	20.11.2012 17:05	MOL File
2821.mol	20.11.2012 17:05	MOL File
2822.mol	20.11.2012 17:05	MOL File
2823.mol	20.11.2012 17:05	MOL File
2824.mol	20.11.2012 17:05	MOL File
2825.mol	20.11.2012 17:05	MOL File
2826.mol	20.11.2012 17:05	MOL File
2827.mol	20.11.2012 17:05	MOL File
2828.mol	20.11.2012 17:05	MOL File
2829.mol	20.11.2012 17:05	MOL File
2830.mol	20.11.2012 17:05	MOL File
list.tsv	20.11.2012 17:05	TSV



list.tsv - Notepad

File Edit Format View Help


l-(3-Chlorophenyl)piperazine	2818.mol
1-[3-(trifluoromethyl)phenyl]piperazine	2819.mol
1-Benzylpiperazine	2820.mol
Amitriptyline	2821.mol
Amphetamine	2822.mol
Benzoyllecgonine	2823.mol
Cocaine	2817.mol
Dextromethorphan	2824.mol
EDDP (2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidin)	2825.mol
Ephedrine	2758.mol
Ketamine	2826.mol
Mephedrone	2827.mol
Methadone	2828.mol
Methamphetamine	2829.mol
Naltrexone	2830.mol

This is the list that connects
structure and record

Final Tasks for MassBank records:

Create “reodata.zip” and “moldata.zip” (e.g. with 7Zip)

- This wasn't automated as this requires additional packages for Windows



These zip files are now ready for upload to MassBank

- If these are new records for NORMAN MassBank
 - Email these to massbank@normandata.eu
- If these are private records, you can upload to your own MassBank
 - I'll give a few screenshots in the next few slides

Uploading MassBank Records

Example with “uchem-massbank” screenshots

Eawag Uchem-MassBank

Spectrum Search



Quick Search



Substructure Search




Browse Page



Peak Search



Spectral Browser



Record Index



Uploading MassBank Records

Example with “uchem-massbank” screenshots



MassBank Administration Tool

Main Menu

[Record Validator](#)

[Instrument Editor](#)

[Record Registration](#)

[Record List](#)

[Structure Registration](#)

[Structure List](#)

[Validator](#)

[SQL File Generator](#)

[File Upload](#)

[Record List Generator](#)

[Version Information](#)

[Database Manager](#)

[HOME](#)

Uploading MassBank Records

Record validation

ADMIN **Record Validator** | Record Registration | Structure Registration | File Upload | Sql File Generator
 MENU | Instrument Editor | Record List | Structure List | Validator | Record List Generator

Record Validator

Database : UchemEZ Record Version : 2 1 (old record version)

info : validation archive is [reodata.zip] or [*.msbk].

0 ok, 210 warn, 0 error / 210 files

Name	Status	Details
EA275801.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275802.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275803.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275804.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275805.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275806.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275807.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275808.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275809.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275810.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.
EA275811.txt	warn	value of required item [CH\$COMPOUND_CLASS:] is not compound class format.

Technically, this should define if the compound is a natural product or not; there is no way to automate this (unambiguously) – hence this warning
 We print “N/A: Environmental Standard

Uploading MassBank Records

Record Registration

uchem-massbank/MassBank/mbadmin/RecordRegist.jsp?db=UchemEZ

ADMIN | **Record Validator** | Record Registration | Structure Registration | File Upload | Sql File Generator | Version Information | HO

MENU | Instrument Editor | Record List | Structure List | Validator | Record List Generator | Database Manager

ADMIN | Record Validator | **Record Registration** | Structure Registration | File Upload | Sql File Generator

MENU | Instrument Editor | Record List | Structure List | Validator | Record List Generator

Record Registration

Database : UchemEZ Record Version : 2 1 (old record version)


Record Archive :

** please specify your [recdata.zip] or [*msbk].*

info : 210 record registered.
info : Done.

Uploading MassBank Records

Structure Registration



uchem-massbank/MassBank/mbadmin/StructureRegist.jsp?db=UchemEZ

ADMIN MENU | Record Validator | Record Registration | **Structure Registration** | File Upload | Sql File Generator

ADMIN MENU | Record Validator | Record Registration | **Structure Registration** | File Upload | Sql File Generator | Validator | Record List Generator

ADMIN MENU | Instrument Editor | Record List | Structure List | Validator | Record List Generator

Database: **Structure Registration**

Database : UchemEZ

Structure Archive : Browse... Registration

** please specify your [moldata.zip] or [*.*sbk].*

info : 15 molfile registered.
info : Done.

Records Registered

Browse MassBank to see what they look like!

Record Index

[mass calculator](#)

[Home](#) | [Spectrum](#) | [Quick](#) | [Peak](#) | [Substructure](#) | [Browser](#) | [Browse](#) | [Index](#) | MassBank ID:

Contributor : [CASMI \(42\)](#) [Chubu Univ. \(2,628\)](#) [EA Uchem Orbi \(5,132\)](#)
[EA Uchem Orbi Test \(641\)](#) [EQ Uchem Q Ex \(1,262\)](#) [Eawag Uchem Adducts \(153\)](#)
[Fukuyama Univ. \(340\)](#) [IMM, CAMS & PUMC, China \(192\)](#) [IPB Halle \(528\)](#)


<input type="checkbox"/>	<input checked="" type="checkbox"/> Amphetamine	14 spectra	C9H13N	135.10480	
					
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 15%; R=15000; [M+H]+				EA282208
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 15%; R=7500; [M+H]+				EA282202
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 30%; R=15000; [M+H]+				EA282209
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 30%; R=7500; [M+H]+				EA282203
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 35%; R=15000; [M+H]+				EA282214
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 35%; R=7500; [M+H]+				EA282201
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 45%; R=15000; [M+H]+				EA282210
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 45%; R=7500; [M+H]+				EA282204
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 60%; R=15000; [M+H]+				EA282211
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 60%; R=7500; [M+H]+				EA282205
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 75%; R=15000; [M+H]+				EA282212
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 75%; R=7500; [M+H]+				EA282206
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 90%; R=15000; [M+H]+				EA282213
<input type="checkbox"/>	LC-ESI-ITFT; MS2; 90%; R=7500; [M+H]+				EA282207

MassBank Record: EA282210

[Home](#) | [Spectrum](#) | [Quick](#) | [Peak](#) | [Substructure](#) | [Browser](#) | [Browse](#) | [Index](#) | MassBank ID:


Amphetamine; LC-ESI-ITFT; MS2; 45%; R=15000; [M+H]⁺

Mass Spectrum



<< < > >> show all m/z

Chemical Structure



ACCESSION: EA282210

RECORD_TITLE: Amphetamine; LC-ESI-ITFT; MS2; 45%; R=15000; [M+H]⁺

DATE: 2012.11.20

AUTHORS: M. Stravs, E. Schymanski, H. Singer, Eawag

LICENSE: CC BY-SA

COPYRIGHT: Copyright (C) Eawag, 2012

COMMENT: CONFIDENCE standard compound

COMMENT: INTERNAL_ID 2822

CH\$NAME: Amphetamine

CH\$NAME: Amfetamine

Take-Home Messages: RMassBank Demo

That was a 60+-slide introduction to RMassBank

- Install (and update) RMassBank and necessary programs
- Run RMassBank for trial data
 - Get a basic understanding for the workflow
- Understand the manual checking required
 - “Fail peaks” – checking these in raw data
 - Automatic annotation with CTS – searching/editing this data
- Generate records with RMassBank
- A quick impression how to upload these records

It will be impossible to remember everything!

- vignette(“RMassBank”) – this contains an extensive explanation
- Stravs et al. 2012, *J. Mass Spectrom.*, DOI: 10.1002/jms.3131

Acknowledgements

- Marie-Curie Post Doctoral Fellowship (E. Schymanski),
- Eawag Discretionary Funds (M. Stravs)
- Coauthors: **Michael Stravs**, Heinz Singer, Juliane Hollender
- R & MassBank help: Steffen Neumann, Michael Gerlich, Erik Müller, (IPB)
- Discussions: Matthias Ruff, Martin Loos (Eawag);
Tobias Schulze, Martin Krauss, Werner Brack (UFZ)
- MassBank & Naming Rights: Prof. Takaaki Nishioka
- NORMAN Association
- And thank you all for listening!



Any Questions?

massbank@normandata.eu (Tobias – record upload)

massbank@eawag.ch (RMassBank: E. Schymanski, M. Stravs)

<http://metabolomics-forum.com/viewforum.php?f=29>

emma.schymanski@eawag.ch

DOI: 10.1002/jms.3131

All Details Contained Within!