

Analytical challenges for the analysis of biocides in aqueous and solid environmental matrices

Arne Wick, Kathrin Broeder, Michael Schluesener, Thomas Ternes
Federal Institute of Hydrology (BfG), Koblenz, Germany




Environmental monitoring of biocides in Europe - from prioritisation to measurements
Workshop, Berlin, 05-06.11.2012


Biocides – main groups





- Disinfectants and general biocidal products
e.g. Human hygiene biocidal products




- Preservatives
e.g. film preservatives



- Pest control, e.g. insecticides




- Other biocidal products, e.g. antifouling



Biocides – Ecotoxicity

Periphyton



Mohr et al. (2008),
Aquatic Toxicol.,
90(2), 109-120

$EC_{50} = 340 \text{ ng/L}$
 $EC_{10} = 10 \text{ ng/L}$


Myriophyllum verticillatum



Mohr et al. (2009),
Environ. Sci. Technol.,
43, 6838-6843

$EC_{50} = 210 \text{ ng/L}$

→ proposed
AA-EQS: 2.5 ng/L



Example: Irgarol

Copepods



Mohr et al. (2008),
Aquatic Toxicol.,
90(2), 109-120

$EC_{50} = 90 \text{ ng/L}$


Oehlmann and Watermann
(2005), UBA Dessau

Radix balthica



$EC_{10} = 32 \text{ ng/L}$

Emission routes of biocides



adapted from Gerecke et al. (2002), 48, 307-315

Analytical requirements

- Broad compound spectrum with different physico-chemical properties
 - positively charged:** *climbazole, imazalil*
 - negatively charged:** *mecoprop, 2,4-D*
 - polar:** *carbendazim* ($\log K_{OW} = 1.5$)
 - non-polar:** *triclocarban* ($\log K_{OW} = 5.1$)
- Complex matrices (surface water, wastewater, sludge, sediment)
- Sensitivity and reproducibility (**e.g. LOQ Irgarol < 2.5 ng/L**)
- High throughput

➔ **Methods based on LC Tandem MS detection**

➔ **Problem: Matrix effects**
→ **often severe ion suppression using ESI**



ESI source

Reduction:

Sample preparation:

- Enrichment volume for SPE
- SPE cartridge
- SPE cartridge wash step
- Clean-up


LC-MS/MS:

- Injection Volume
- Chromatographic conditions
- Flow rate to MS (split)
- **Ionization source (APCI vs. ESI)**


Compensation:

- **Use of labeled surrogate standards (isotope dilution technique)**
- Use of standard addition

Analytical method – water samples



Analytical method – solid samples



Pressurized Liquid Extraction (PLE)



ASE Parameter: 80°C/100bar
Different extraction solvents

Measurement with LC-MS/MS

~ 45 biocides and pesticides, + benzothiazoles and polar UV-filters

- Isothiazolones (e.g. Benzisothiazolone, BIT, Octylisothiazolone, OIT)
- Imidazole fungicides (e.g. climbazole, imazalil)
- Triazole fungicides (e.g. propiconazole, tebuconazole)
- Carbamate fungicides (e.g. carbendazim, IPBC)
- Phenyl urea herbicides (e.g. diuron, isoproturon)
- Triazines (e.g. terbutryn, irgarol)
- Bacteriocides (e.g. triclosan, triclocarban)
- Insecticides (e.g. DEET, imidacloprid)


Matrix: groundwater, surface water, wastewater, partly for sludge, sediments and soils

LOQ_{surface water}: 0.5 – 50 ng/L, Accuracy: 80-120%, Precision < 25%

Determined for every analytical run and matrix!!

Determination of matrix effects

Post-extraction spike




Matrix effect (ME)


= Absolute recovery of a post-extraction spike

$$\text{ME(\%)} = \frac{\text{peak area}_{\text{post-extraction spike}} - \text{peak area}_{\text{blank}}}{\text{peak area}_{\text{external standard}}} \cdot 100$$







Matrix effects (raw wastewater)



Relative Recovery



Example: Thiabendazole



S/N = 35, ~ 4 ng/L

→ LOQ ~ 1 ng/L

S/N = 9, ~ 4 ng/L


→ LOQ ~ 5 ng/L

Direct injection

No enrichment prior to measurement with an API 5500 Tandem MS

→ No losses during SPE, High Throughput

→ Less matrix effects !



→ Possible drawbacks: Higher LOQs, higher blank values


Extraction of solid samples

Many biocides are predominantly sorbed:

Quarternary ammonium compounds (QACs), triclosan, triclocarban, chlorophene, many conazoles (e.g. climbazole, ketoconazole, miconazole), etc.

→ Analysis in particulate matter, sediments, biota, soils and sludge is required

→ Challenge: Extractability depends on the physico-chemical properties of the compounds (e.g. charge, K_{OW} , etc.) and the properties of the solid matrix (cation exchange capacity, pH, TOC, etc.) → *no certified reference materials*




Extraction of solid samples

Recovery [%] of climbazole using self-prepared reference materials

	MeOH/H ₂ O	MeOH	MeOH/Acetone
<i>soil</i>	79	60	32
<i>sediment</i>	78	51	8
<i>sludge</i>	50	71	20


Recovery using freeze-dried soil spiked with climbazole

directly prior to extraction with MeOH/acetone: ~ 100 %




Concentration variability

Concentration [ng/L]



Terbutryn (n = 11)

AA EQS = 65 ng/L



Irgarol (n = 11)

M1 in the same conc. range

AA EQS = 2.5 ng/L

Conclusions

Many biocides can be determined simultaneously in various environmental matrices down to the low ng/L and ng/g level. Ionization by **APCI** can be a valuable alternative, especially if labeled **surrogate standards** are missing

*

A compensation of **matrix effects** is the prerequisite for analytical accuracy and reproducibility.

*

The **relative recoveries** and the **LOQs** should be determined for each sample series and each sample matrix.

*

Extraction methods for solid matrices should be validated for each individual matrix (e.g. soils, sediments, sludge, biota). Self-prepared reference materials are an option for testing extractability.

*

A sufficient **time resolution for sampling** is crucial for assessing annual average and maximum concentrations of biocides

An aerial photograph of a wide river flowing through a valley. Several barges are visible on the water, including one with a red and blue hull and another with a black deck. The riverbanks are lined with dense green trees and a road. In the background, a town is nestled in a valley, and rolling green hills rise under a clear blue sky.

Thank you for your attention

Questions ?